

# Product datasheet BIOCRYL® C

## 1 Description of the product / intended use / company

Trade name:	<b>BIOCRYL® C</b>
Intended use:	Thermoplastic foils for the manufacture of appliances and splints for dental and orthodontic treatments.
Manufacturer / Supplier:	Scheu Dental GmbH Am Burgberg 20 58642 Iserlohn Germany Tel. 0049 2374 9288-0 <a href="mailto:qmra@scheu-dental.com">qmra@scheu-dental.com</a>

## 2 Composition / information on ingredients

### Chemical characterization

CAS-Number:	9011-14-7
Designation:	Polymethylmethacrylat (PMMA)

## 3 Handling and storage

### Handling

Information for safe handling:	When using do not eat, drink or smoke. Provide suction extractors if dust is formed.
--------------------------------	---

### Storage

Storage conditions:	Store dry and dark at 5 - 35 °C
---------------------	---------------------------------

## 4 Physical, chemical, mechanical and biological properties

### 4.1 General properties

Properties	Guideline	Value
Form	-	Solid
Colour	-	Depending on the version, transparent or colored
Odour	-	Odourless
Density	ISO 1183	1,19 g/cm <sup>3</sup>
Water absorption after 24 h at 23°C	-	-

# Product datasheet BIOCRYL® C

## 4.2 Mechanical properties

Properties	Guideline	Value
Tensile strength	ISO 527	76 MPa
Tensile modulus	ISO 527	3200 MPa
Elongation at break	ISO 527	5,5 %
Charpy Impact Strength @ 23 °C	ISO 179/1eU	20 kJ/m <sup>2</sup>

## 4.3 Thermal properties

Properties	Guideline	Value
Glass transition temperature	ISO 11357	112 °C
Vicat softening point	ISO 306	103 °C

## 4.4 Biological properties / Biocompatibility

A biological evaluation was carried out for the material in accordance with DIN EN ISO 10993-1. It meets the requirements regarding biocompatibility for medical devices.

## 5 Stability and reactivity

### Thermal decomposition / Conditions to avoid:

No decomposition by intended use.

## 6 Disposal

The raw material can be recycled after separation or disposed of like commercial or household waste.

The above information is provided to the best of our knowledge. In addition to the information provided, it is recommended to check the suitability of the products individually regarding the respective application and associated procedures and purposes.

Serious incidents in connection with the use of the product shall be reported immediately to the manufacturer ([gmra@scheu-dental.com](mailto:gmra@scheu-dental.com)) and the competent authority.