

Product datasheet DURASOFT® seal

1 Description of the product / intended use / company

Trade name:	DURASOFT® seal
Intended use:	Thermoplastic foils for the manufacture of appliances and splints for dental and orthodontic treatments.
Manufacturer / Supplier:	Scheu Dental GmbH Am Burgberg 20 58642 Iserlohn Germany Tel. 0049 2374 9288-0 qmra@scheu-dental.com

2 Composition / information on ingredients

Chemical characterization

CAS-Number:	75701-44-9
Designation:	Thermoplastic Polyurethane (TPU)

3 Handling and storage

Handling

Information for safe handling:	When using do not eat, drink or smoke. Provide suction extractors if dust is formed.
--------------------------------	---

Storage

Storage conditions:	Store dry and dark at 5 - 35 °C
---------------------	---------------------------------

4 Physical, chemical, mechanical and biological properties

4.1 General properties

Properties	Guideline	Value
Form	-	Solid
Colour	-	Depending on the version, transparent, blue or magenta
Odour	-	Odourless
Density	ASTM D 792	1,20 g/cm ³
Water absorption after 24 h at 23°C	ASTM D 570	0,19 %

Product datasheet DURASOFT® seal

4.2 Mechanical properties

Properties	Guideline	Value
Tensile strength	ASTM D 638	69 MPa
Tensile modulus	ASTM D 638	2100 MPa
Elongation at break	ASTM D 638	140 %
Rockwell hardness (R scale)	ASTM D 758	123

4.3 Thermal properties

Properties	Guideline	Value
Vicat softening point	ASTM D 1525	109 °C

4.4 Biological properties / Biocompatibility

A biological evaluation was carried out for the material in accordance with DIN EN ISO 10993-1. It meets the requirements regarding biocompatibility for medical devices.

5 Stability and reactivity

Thermal decomposition / Conditions to avoid:

No decomposition by intended use.

6 Disposal

The raw material can be recycled after separation or disposed of like commercial or household waste.

The above information is provided to the best of our knowledge. In addition to the information provided, it is recommended to check the suitability of the products individually regarding the respective application and associated procedures and purposes.

Serious incidents in connection with the use of the product shall be reported immediately to the manufacturer (gmra@scheu-dental.com) and the competent authority.